


EasyCell[®] assistant

Remote Review Software



Add workstations for greater efficiency and collaboration

Now technologists can review cell images from EasyCell assistants from anywhere. This allows for increased throughput for reviews and the added convenience of location independence.

Greater staffing efficiency

Assign multiple technologists to review EasyCell cell images and sign out patient samples. EasyCell Remote permits easy addition of staff for review of samples when needed.

Access to expert opinions

Share images with colleagues located at distant centers of excellence. Images from difficult cases can be conveniently sent to second opinions.


Monitor quality in the lab

View patient results from all areas of the lab. Supervisors and administrators can review sample results and assure consistency in reporting.

Instant clinician access to image results

Decisions regarding patient treatment may be accelerated with rapid viewing of images from any pathologist's office. Access to EasyCell image files permits comparison of results over the history of the case.

EasyCell Remote Software



Key Features

Centralized Expertise

- Run EasyCell® Remote Review Software alongside other applications on any existing PC
- Review preclassifications from one or more EasyCell assistants
- Access all patient results
- Display associated patient information (age, gender, etc.)
- Up to 16 concurrent viewers

Technical Specifications

Recommended PC specifications for EasyCell Remote

- Windows XP, Windows Vista, or Windows 7
- 64 MB video card
- 512 MB RAM
- 50 MB free disk space
- Pentium 4 or higher

MEDICA