



EasyCell[®] assistant

Competency Review Software

Evaluate technologist's expertise in classifying white cells and red cell morphology

- Standardize cell differentials
- Enhance quality of results
- Achieve confidence in cell reporting

Built-in competency review

No additional software or internet connection required.
(optional activation available)

Easy Integration

Instructors can add new, reviewed, or reported slide files to the competency review set from their laboratory's files. Competency review is easily integrated into the daily routine.

MEDICA

Competency Review Software

An instructor selects specific slide files and groups them into a test set. After creating an answer key, the instructor gives the test to one or more students for training purposes or for periodic recertification.

The software intentionally provides confidentiality with proficiency evaluation. Individual skill levels are recognized.

The instructor generates reports showing the student's performance on the test.

Both instructor and student are able to review incorrect answers.

Student Name: Notartomaso, Claire

Instructor Name: Hematology Supervisor

Report Date: 03/27/2013

WBC Classifications

Categories	Instructor Count	Student Count	Matches
Seg	75	74	74
Band			
Lymph	23	20	20
Mono	9	11	9
Eo	3	2	2
Baso	1	1	1
Variant Lymph	5	8	4
Meta			
Myelo			
Promyelo			
Blast			
Plasma			
Other			
NRBC			
Smudge	1	1	0
Unknown	10	9	8

WBC Morphology	Instructor Count	Student Count
Auer Rods	0	0
Döhle Bodies	0	0
Hypersegmentation	1	0
Toxic Granulation	0	0
Vacuolization	0	0

RBC Morphology

Categories	Instructor Count	Student Count
Anisocytosis	0	0
Macrocytes	0	0
Microcytes	0	0
Poikilocytosis	0	1
Acanthocyte	0	0
Burr	1	0
Howell-Jolly Bodies	0	0
Ovalocyte	1	0
Pappenheimer Bodies	0	0
Schistocyte	1	1
Sickle	0	0
Spherocyte	0	0
Target Cell	0	0
Hypochromic	0	0
Polychromatic	0	0
Basophilic Stippling	0	0
Rouleaux	1	0
Overall Normal	0	0

Platelet Results

Platelet Estimate	Instructor Count	Student Count
Normal	1	0
Increased	0	0
Decreased	0	0

Platelet Morphology	Instructor Count	Student Count
Clumped	0	0
Giant	0	0

Additional Comments :

Great Job Claire!

